




Bair Hugger® Therapy Cardiothoracic Blanket Options

Models 630, 635, 637, 645

Designed to warm patients during surgery so they arrive to the ICU at normothermic temperatures




Sterile Cardiac Blanket

- Warms the lower body during surgery to reduce core to peripheral heat loss
- Positions easily in the sterile field
- Tapes securely to the waistline
- Clear groin window allows access to femoral artery (5¾ in. x 15¾ in.) (15 x 40 cm)
- Foot drape minimizes risk of thermal injury to feet and lower legs
- 36 inch sterile hose allows the placement of the temperature management unit away from the sterile field



Model 630



Full Access & Sterile Full Access Underbody Blankets

- Pass-through slit design allows for flexible patient positioning and drawsheet use
- Fluid outlets minimize pooling of fluids on the surface of the blanket
- Soft, radiolucent materials accommodate imaging requirements
- Perforated area at head of blanket can be removed to accommodate a full prone position
- Adhesive strips and tuck flaps secure the blanket to the O.R. / procedure table
- Resealable hose ports at either end of the blanket provide options for hose placement
- One clear plastic drape (included) helps retain warm air around the patient's head




Model 635 (non-sterile)/637 (sterile)



Cardiothoracic Blanket Options

Models 630, 635, 637, 645




Sterile Cardiac Access Blanket

- Warms the patient's periphery during surgery
- Radiolucent materials will not interfere with imaging requirements
- Tapes securely to skin so the blanket can provide localized warming throughout the surgical procedure
- Allows full access to both legs
- Insulated foot pouch minimizes risk of thermal injury to feet



Model 645

Product Offerings



Specifications:

Sterile Cardiac Blanket

Model 630
 Size: 60" x 36" (152 x 91 cm)
 Weight: 8.0 oz (227 g)
 Part# - 63000
 Units/Case - 5

Full Access Underbody Blanket

Model 635
 Size: 84" x 36" (221 x 91 cm)
 Weight: 7.0 oz (198 g)
 Drape: 24" x 48" (61 x 122 cm)
 Part# - 63500
 Units/Case - 5

Sterile Cardiac Access Blanket

Model 645
 Size: 9" x 108" (23 x 274 cm)
 Weight: 4.0 oz (113 g)
 Part# - 64500
 Units/Case - 5

Sterile Full Access Underbody Blanket

Model 637
 Size: 84" x 36" (221 x 91 cm)
 Weight: 7.0 oz (198 g)
 Drape: 24" x 48" (61 x 122 cm)
 Part# - 63700
 Units/Case - 5



©2008 Arizant Healthcare Inc. All rights reserved.
 10393 West 70th Street, Eden Prairie, MN 55344 USA
 tel: 952-947-1200 | 800-733-7775
 fax: 952-947-1400 | 800-775-0002
 web: arizanthealthcare.com | bairhugger.com

Bair Hugger, Arizant, Arizant Healthcare, and the Bair Hugger and Arizant logos are trademarks of Arizant Healthcare Inc., registered in the U.S. Patent & Trademark Office and in other countries.

601214D 11/08

To order, please contact your local Arizant Healthcare® sales representative or call Customer Service at 1-800-733-7775.

