



Bair Hugger[®] Therapy Temperature Management Unit

Model 750



Performance. Precision. Ease of Use.

Model 750



The Next Generation of Temperature Management Technology



Performance

- High airflow for greater heat transfer to the patient
- Hose end temperature sensing ensures consistent temperatures

Precision

- Microprocessor technology regulates temperature settings within $\pm 1.5^{\circ}\text{C}$ for accurate temperature delivery
- Three independent sensors continuously monitor system performance

Ease of Use

- Advanced user interface provides finger-tip control of temperature settings for one-touch temperature selection
- No need to open the unit for over-temperature tests or calibration

Specifications

Dimensions

12"h x 10"w x 13.5"d (30h x 25w x 34d cm)

Weight

15.5 lb (7 kg)

Operating Temperatures

High: $43^{\circ} \pm 1.5^{\circ}\text{C}$ ($109.4^{\circ} \pm 2.7^{\circ}\text{F}$)
Med: $38^{\circ} \pm 1.5^{\circ}\text{C}$ ($100.4^{\circ} \pm 2.7^{\circ}\text{F}$)
Low: $32^{\circ} \pm 1.5^{\circ}\text{C}$ ($89.6^{\circ} \pm 2.7^{\circ}\text{F}$)

Leakage Current

Meets regulatory standards for leakage current

Filter

High-efficiency 0.2 μm filter

Device Ratings

110-120 VAC, 50/60 Hz, 11.7 Amperes
220-240 VAC, 50/60 Hz, 7.2 Amperes
100 VAC, 50/60 Hz, 14.5 Amperes

Ordering Information

Bair Hugger Model 750 Temperature Management Unit

To order, please contact your Arizant Healthcare Inc. representative or call 1-800-733-7775.


**Arizant
Healthcare**

©2008 Arizant Healthcare Inc. All rights reserved.
10393 West 70th Street, Eden Prairie, MN 55344 USA
tel 952-947-1200 | 800-733-7775
fax 952-947-1400 | 800-775-0002
web arizanthealthcare.com | bairhugger.com

Bair Hugger, Arizant and the Bair Hugger and Arizant logos are trademarks of Arizant Healthcare Inc. registered in the U.S. Patent & Trademark Office.