

BAIR HUGGER THERAPY & RANGER SYSTEM TOGETHER

Complete Pediatric Temperature Management

RANGER SYSTEM PEDIATRIC FLUID WARMING



Warming Unit - Model 24500
Set point temperature: 41°C
Flow rates: KVO – 30,000 mL/hr
Over temperature set points:
Primary: 43°C
Secondary: 46°C
Power: 110VAC/220-240VAC
Weight: 7 lbs. 7 oz.
Dimensions: 7.5" w x 4.5" h x 10" d
 (19 x 11 x 25 cm)



Pediatric/Neonate Disposable Set Model 24400
Flow rates: KVO - 6,000 mL/hr
Priming volume: 20 mL
Components: Heat exchanger, bubble trap with air aspiration port, patient line
Pressure infusion: Up to 300 mmHg
Sterile: EtO
Cassette: DEHP free cassette
Use: Single use only
Units/case: 10



Pediatric/Neonate Disposable Set with Fluid Aspiration Port Model 24450
Flow rates: KVO - 6,000 mL/hr
Priming volume: 20 mL
Components: Heat exchanger, bubble trap with air aspiration port, fluid aspiration port, patient line
Pressure infusion: Up to 300 mmHg
Sterile: EtO
Cassette: DEHP free cassette
Use: Single use only
Units/case: 10

BAIR HUGGER THERAPY PEDIATRIC BLANKETS



Pediatric Underbody Blanket Model 55501
Size: 36 x 33 in (91 x 84 cm)
Weight: 3.0 oz (85 g)
Drape (2): 24 x 24 in (61 x 61 cm)
Units/case: 10



Large Pediatric Underbody Blanket Model 55000
Size: 60 x 32 in (152 x 81 cm)
Weight: 4.8 oz (136 g)
Drape (2): 24 x 24 in (61 x 61 cm)
Units/case: 10



Small Lower Body Blanket Model 53700
Size: 35 x 24 in (89 x 61 cm)
Weight: 1.5 oz (43 g)
Units/case: 10



Pediatric Blanket Model 31000
Size: 60 x 36 in (152 x 91 cm)
Weight: 4.8 oz (134 g)
Units/case: 10

To arrange an evaluation of any of these products, contact your local Arizant Healthcare representative or call 1-800-733-7775.



Ranger
 Blood/Fluid Warming System

COMPLETE PEDIATRIC TEMPERATURE MANAGEMENT

We're not hiding...



Ranger
 Blood/Fluid Warming System



INTRODUCING COMPLETE PEDIATRIC WARMING SOLUTIONS



RANGER* SYSTEM PEDIATRIC FLUID WARMING

Ranger Pediatric/Neonate Blood and Fluid Warming Disposable Sets



Pediatric/Neonate Disposable Set Model 24400



Pediatric/Neonate Disposable Set with Fluid Aspiration Port Model 24450

The unique design of the pediatric/neonate disposable heat exchanger set allows the Ranger system to fully address pediatric fluid warming needs. The sets are designed for low flow, maintenance fluid delivery and aliquot fluid administration where measured fluid volumes are infused via syringe.

BAIR HUGGER* THERAPY PEDIATRIC BLANKETS

Underbody, lower body and full blankets sized just right for your pediatric patients, from neonates to young adults.

Pediatric underbody blankets allow full access to the patient and include two drapes that help retain the warm air that surrounds the intubated patient.



Model 55501



Model 31000



Model 55000



Model 53700

Only underbody blankets feature patented fluid outlets that minimize pooling of fluids on the blanket's surface.

What you're seeking!



Complete Pediatric Temperature Management

You know the importance of keeping your littlest patients warm, so why trust them to anyone but the leader in patient warming? Arizant Healthcare Inc., the creator of forced-air warming, brings you the latest advances in surgical warming for pediatric patients.

Pediatric warming solutions that are designed specifically for pediatric applications—no need to adapt products and equipment for pediatric use.